



Comfort Plus™ Premium Saliva Ejectors

INSTRUCTIONS FOR USE

1. Attach ejector head onto the end of the low volume evacuation tube.
2. Place in mouth and bend to desired position.
3. After use, dispose of the ejector.



Ref.#'s ZWWCP, ZCWCP,
ZWPCPBG, ZWGCPM

 Crosstex International, Inc.
6789 W. Henrietta Road, Rush, NY 14543 USA
585.359.0130 | 631.582.6777 | crosstex.com

 MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany


0050



Not made with natural rubber latex.

Made in USA
ZWCP-IFU-EN / REV B - 03/2022

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.
©2022 Crosstex International, Inc. All rights reserved.

CROSSTEX