




# HU-FRIEDY KEY OPINION LEADER

## ACADEMY OF RICHARD V. TUCKER STUDY CLUBS KIT LIST

### The Academy of Richard V Tucker Study clubs

consist of 75 independently operating, clinical study clubs throughout the United States, Canada, Germany, Italy, and Japan. Each study club consists of approximately 12 members and a mentor. The mentor ensures that the study club session adheres to the Tucker technique. Each session is divided into a clinical component and a critique.



PART CODE	DESCRIPTION	QUANTITY
<b>MT232TM9</b>	232 TruBal Margin Trimmer #9Hdl	
<b>MT233TM9</b>	233 TruBal Margin Trimmer #9Hdl	
<b>MT232TBM9</b>	232 Tucker Margin Trimmer #9Hdl	
<b>MT233TBM9</b>	233 Tucker Margin Trimmer #9Hdl	
<b>CP14/14-09</b>	#14/14-0 DE Hatchet (15-10-16) #9Hdl	
<b>CP15/15-09</b>	#15/15-0 DE Hatchet (15-10-16) #9Hdl	
<b>CP44S9</b>	#44S Off-Angle Cavity Prep #9Hdl	
<b>CP45S9</b>	#45S Off-Angle Cavity Prep #9Hdl	
<b>IM5109</b>	10 Inst. Sig Cassette Green	
<b>KIT-TUCKER</b>	<b>Complete Tucker Kit</b>	

Check One: Quote ☐ Order ☐

Contact 1-800-Hu-Friedy for more information or fax order to 1-800-729-1299.

Name \_\_\_\_\_

Address \_\_\_\_\_

City \_\_\_\_\_

State \_\_\_\_\_

Zip Code \_\_\_\_\_

Phone \_\_\_\_\_

Preferred Dealer \_\_\_\_\_

Dealer Rep Name \_\_\_\_\_

Please check the box if you would like to receive confirmation: ☐



# THE ACADEMY OF RICHARD V TUCKER STUDY CLUBS KIT LIST

**A** MT232TM9

**B** MT233TM9

**C** MT232TBM9

**D** MT233TBM9


**E** CP14/14-09

**F** CP15/15-09

**G** CP44S9

**H** CP45S9

**I** IM5109



## THE ACADEMY OF RICHARD V TUCKER STUDY CLUBS KIT LIST



- A MT232TM9**  
**232 TruBal Margin Trimmer #9Hdl**  
Margin Trimmer used to produce proper bevel on enamel margins. Similar to a hatchet except the blade is curved and the cutting edge angled.



- B MT233TM9**  
**233 TruBal Margin Trimmer #9Hdl**  
Margin Trimmer used to produce proper bevel on enamel margins. Similar to a hatchet except the blade is curved and the cutting edge angled.



- C MT232TBM9**  
**232 Tucker Margin Trimmer #9Hdl**  
Margin Trimmer used to produce proper bevel on enamel margins. Similar to a hatchet except the blade is curved and the cutting edge angled. Margin Trimmer used to produce proper bevel on enamel margins. Similar to a hatchet except the blade is curved and the cutting edge angled.



- D MT233TBM9**  
**233 Tucker Margin Trimmer #9Hdl**  
Margin Trimmer used to produce proper bevel on enamel margins. Similar to a hatchet except the blade is curved and the cutting edge angled.



- E CP14/14-09**  
**#14/14-0 DE Hatchet (15-10-16) #9Hdl**  
Used for cavity preparation: retentive areas, internal line angles and removing hard caries. 1.5mm/1.5mm [15-10-16]



- F CP15/15-09**  
**#15/15-0 DE Hatchet (15-10-16) #9Hdl**  
Used for cavity preparation: retentive areas, internal line angles and removing hard caries. 1.5mm/1.5mm [15-10-16]



- G CP44S9**  
**#44S OFF-ANGLE CAVITY PREP #9HDL**  
Used for cavity preparation: retentive areas, internal line angles and removing hard caries.



- H CP45S9**  
**#45S Off-Angle Cavity Prep #9Hdl**  
Used for cavity preparation: retentive areas, internal line angles and removing hard caries.



- I IM5109**  
**10 Inst. Sig Cassette Green**  
IMS® cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.