

Technical Data Sheet

Crosstex Sparkle® Propyl Paste with Xylitol



<u>PRODUCT</u>	Crosstex Sparkle® Propphy Paste with Xylitol																																
<u>PRODUCT OVERVIEW</u>	Sparkle Prophylaxis Paste is a blend of polishing and cleaning agents, which makes it effective in removing stains and smoothing the surfaces of teeth. Variety of available grits allows you to choose the least abrasive polishing agent to remove biofilm and stain.																																
<u>MATERIALS</u>	<ul style="list-style-type: none">• Glycerin• Water	<ul style="list-style-type: none">• Pumice• Sodium Fluoride	<ul style="list-style-type: none">• FD&C Coloring and flavoring• Xylitol																														
<u>DETAILS</u>	<ul style="list-style-type: none">• Contains 1.23% Fluoride ion• Formulated to maximize stain removal with minimal enamel loss• Advanced formulation provides consistent, splatter-free results• Completely rinses out, leaves no residual grit• Available in multiple grits and flavors																																
<u>ALLERGEN INFORMATION</u>	Product is free of nuts, soy, milk, eggs, gluten, aspartame and corn/corn products and is manufactured in area in which has no exposure to the following items: wheat, eggs, peanut, milk, fish, nut, soybean and shellfish																																
<u>PRODUCT/ PACKAGING</u>	<div>Sparkle® Propphy Paste is available in the following configurations</div> <table><tr><th>Product Code</th><th>Description</th></tr><tr><td>UPC12</td><td>Course grit, Mint, 12 oz/340 g jar</td></tr><tr><td>UPCB</td><td>Course grit, Berrylicious Assortment*, 200/Box</td></tr><tr><td>UPCBG</td><td>Course grit, Bubble Gum, 200/Box</td></tr><tr><td>UPCCH</td><td>Course grit, Cherry, 200/Box</td></tr><tr><td>UPCM</td><td>Course grit, Mint, 200/Box</td></tr><tr><td>UPCOV</td><td>Course grit, Orange Vanilla, 200/Box</td></tr><tr><td>UPCA</td><td>Course grit, Assorted Flavors, 200/Box</td></tr><tr><td>UPMB</td><td>Medium grit, Berrylicious Assortment*, 200/Box</td></tr><tr><td>UPMBG</td><td>Medium grit, Bubble Gum, 200/Box</td></tr><tr><td>UPMCH</td><td>Medium grit, Cherry, 200/Box</td></tr><tr><td>UPMM</td><td>Medium grit, Mint, 200/Box</td></tr><tr><td>UPMOV</td><td>Medium grit, Orange Vanilla, 200/Box</td></tr><tr><td>UPMA</td><td>Medium grit, Assorted Flavors, 200/Box</td></tr><tr><td>UPFBG</td><td>Fine grit, Bubble Gum, 200/Box</td></tr></table>			Product Code	Description	UPC12	Course grit, Mint, 12 oz/340 g jar	UPCB	Course grit, Berrylicious Assortment*, 200/Box	UPCBG	Course grit, Bubble Gum, 200/Box	UPCCH	Course grit, Cherry, 200/Box	UPCM	Course grit, Mint, 200/Box	UPCOV	Course grit, Orange Vanilla, 200/Box	UPCA	Course grit, Assorted Flavors, 200/Box	UPMB	Medium grit, Berrylicious Assortment*, 200/Box	UPMBG	Medium grit, Bubble Gum, 200/Box	UPMCH	Medium grit, Cherry, 200/Box	UPMM	Medium grit, Mint, 200/Box	UPMOV	Medium grit, Orange Vanilla, 200/Box	UPMA	Medium grit, Assorted Flavors, 200/Box	UPFBG	Fine grit, Bubble Gum, 200/Box
Product Code	Description																																
UPC12	Course grit, Mint, 12 oz/340 g jar																																
UPCB	Course grit, Berrylicious Assortment*, 200/Box																																
UPCBG	Course grit, Bubble Gum, 200/Box																																
UPCCH	Course grit, Cherry, 200/Box																																
UPCM	Course grit, Mint, 200/Box																																
UPCOV	Course grit, Orange Vanilla, 200/Box																																
UPCA	Course grit, Assorted Flavors, 200/Box																																
UPMB	Medium grit, Berrylicious Assortment*, 200/Box																																
UPMBG	Medium grit, Bubble Gum, 200/Box																																
UPMCH	Medium grit, Cherry, 200/Box																																
UPMM	Medium grit, Mint, 200/Box																																
UPMOV	Medium grit, Orange Vanilla, 200/Box																																
UPMA	Medium grit, Assorted Flavors, 200/Box																																
UPFBG	Fine grit, Bubble Gum, 200/Box																																
	*Includes: Raspberry, Blueberry, Strawberry																																
<u>STORAGE REQUIREMENTS</u>	n/a																																
<u>SHELF LIFE</u>	2 years																																
<u>REGULATORY INFORMATION</u>	Product Code: EJR Classification: FDA Class 1 Regulation: Medical Device Listing: E646228																																

*Includes: Raspberry, Blueberry, Strawberry