

Technical Data Sheet

Crosstex Sparkle® FREE Propyl Paste with Xylitol



PRODUCT

Crosstex Sparkle® FREE Propyl Paste with Xylitol

PRODUCT OVERVIEW

Sparkle Free Prophylaxis Paste is a blend of polishing and cleaning agents, which makes it effective in removing stains and smoothing the surfaces of teeth. Variety of available grits allows you to choose the least abrasive polishing agent to remove biofilm and stain.

MATERIALS

- Glycerin
- Water
- Pumice
- Xylitol

DETAILS

- Prophy Paste in a unique fluoride/dye free formula
- Formulated to maximize stain removal with minimal enamel loss
- Advanced formulation provides consistent, splatter-free results
- Completely rinses out, leaves no residual grit
- Available in multiple grits and flavors

ALLERGEN INFORMATION

- Product is free of nuts, soy, milk, eggs, gluten, aspartame and corn/corn products and is manufactured in area in which has no exposure to the following items: wheat, eggs, peanut, milk, fish, nut, soybean and shellfish

PRODUCT/ PACKAGING

Sparkle® Free Propyl Paste is available in the following configurations

Product Code	Description
UPS FCC	Course grit, Cinnamon, 200/Box
UPS FCF	Course grit, Fruity, 200/Box
UPS FCS	Course grit, Spearmint, 200/Box
UPS FCWC	Course grit, White Chocolate, 200/Box
UPS FMC	Medium grit, Cinnamon, 200/Box
UPS FMF	Medium grit, Fruity, 200/Box
UPS FMS	Medium grit, Spearmint, 200/Box
UPS FMWC	Medium grit, White Chocolate, 200/Box

STORAGE REQUIREMENTS

n/a

SHELF LIFE

2 years

REGULATORY INFORMATION

Product Code: EJR
Classification: FDA Class 1
Regulation:
Medical Device Listing: E646228