


Crosstex™ Non-Woven Headrest Covers

REF LOCB / LOCL / LOC / L3CB / L3C

INSTRUCTIONS FOR USE

- 1. Use an approved surface disinfectant as recommended by the chair manufacturer to clean the headrest.
- 2. Hold cover with opening on the bottom.
- 3. Place cover over headrest until headrest is completely covered.



Crosstex International, Inc., a Cantel Medical Company
 6789 W. Henrietta Road, Rush, NY 14543 USA
 585.359.0130 | 631.582.6777 | crosstex.com



MDSS GmbH
 Schiffgraben 41
 30175 Hannover, Germany



Not made with natural rubber latex.

Made in China
 NWHEADCVR-IFU-EN / REV C - 05/2021

All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.
 ©2021 Crosstex International, Inc. All rights reserved.