

Cup Dispenser

REF PCC

INSTRUCTIONS FOR USE

1. Pull off top of the dispenser.
2. Load cups with opening pointing up. Dispenser holds 50 cups.
3. Replace top.



CROSSTEX



Crosstex International, Inc.
6789 W. Henrietta Road, Rush, NY 14543 USA
585.359.0130 | 631.582.6777 | crosstex.com

EC REP

MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



MD

Made in Sweden

PCC-IFU-EN / REV B - 06/2021

All product names are trademarks of Crosstex International, Inc., its affiliates or related companies, unless otherwise noted.
©2021 Crosstex International, Inc. All rights reserved.

Not made with natural rubber latex.