

Surgical Aspirator Tips

REF ZSABL / ZSAWH / ZSAGR
INSTRUCTIONS FOR USE

Surgical Aspirator Tips with small, standard or large openings and color-coded according to size. They can fit 11 mm suction holders.

Indication: Blood and fluid suction during surgeries.

1. Connect the Surgical Aspirator Tip to the high volume evacuation hose.
2. Discard aspirator after use.



Small Tip (.0625 in / .16 cm) Blue
Standard Tip (.125 in / .32 cm) White
Large Tip (.25in / .64 cm) Green

CROSSTEX



Manufactured for: Crosstex International, Inc.
6789 W. Henrietta Road, Rush, NY 14543 USA
585.359.0130 | 631.582.6777 | HuFriedyGroup.com



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



0050



Not made with natural rubber latex.

Made in China

ZSA-IFU-EN / REV B - 04/2022

In case of any serious incident that has occurred and involves our devices within the European Union, please report it to Hu-Friedy Mfg. Co., LLC and your national competent authority of the member state.

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.

©2022 Hu-Friedy Mfg. Co., LLC. All rights reserved.