

# T-Handle Light Sleeves

REF BCLST

## INSTRUCTIONS FOR USE

- 1. Locate the arrow on the plastic sleeve. This arrow points to the die-cut.
- 2. Slide plastic sleeve over the handle.



Crosstex International, Inc., a Cantel Medical Company  
6789 W. Henrietta Road, Rush, NY 14543 USA  
585.359.0130 | 631.582.6777 | crosstex.com



MDSS GmbH  
Schiffgraben 41  
30175 Hannover, Germany



Not made with natural rubber latex.

Made in China

BCLST-IFU-EN / REV B - 06/2021

All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.

©2021 Crosstex International, Inc. All rights reserved.