

# Dental Waterline Testing



**The steps you take for infection prevention are important, and that extends to the safety of your dental unit waterlines.**

The narrow tubing can be a breeding ground for a variety of dangerous bacteria, and the consequences of inaction can be dire for your patients and your practice. But treatment alone can't provide you the assurance you need. To be confident, you need to know your water is safe and compliant.

HuFriedyGroup offers both mail-in and in-office waterline testing options. Use these test methods in combination as part of a holistic protocol for waterline maintenance.



## COMPREHENSIVE SUPPORT

**We provide knowledgeable, professional resources and guidance to help you understand your testing results and support your overall water program.**

## 24-Hour In-Office Dental Waterline Test Kit

### FEATURES AND BENEFITS

- **Quick Results**—Detects microbial levels that exceed 500 CFU/mL within 24 hours at room temperature.
- **Simple Process**—Collect water samples directly into provided vials, then just seal and shake.
- **Immediate Neutralization**—Agent inside the vial stops action of any antimicrobial used to treat water, ensuring a more accurate result.
- **Easy-to-Read**—Determine pass/fail with clearly identifiable color change results.



### In-Office

24-Hour In-Office Dental Waterline Test Kit - 4 Vials

**IMS-1004**

The 24-Hour In-Office Dental Waterline Test Kit is for use as a screening tool to provide a pass/fail indication based upon the current EPA standard of  $\leq 500$  CFU/mL. If exact bacteria counts are needed, we recommend testing with our Mail-in Dental Waterline Test Service. Reach out to HuFriedyGroup water support with questions.

## Dental Waterline Test Service

### FEATURES AND BENEFITS

- **Third-Party ISO-Certified Lab**—Testing is conducted by an ISO 17025:2017 accredited laboratory with unbiased results for external validation.
- **Exclusive Platform**—Online member portal provides a single location to track your water samples; retrieve and historically track test results.
- **Simple Process**—Everything that you need – test vials, shipping materials and pre-paid shipping – is included with the kit. Simply collect water samples from your facility and send to the lab.



### Mail-In

Dental Waterline Test Kit, HPC Lab Test - 4 vials

**IMS-2004**

Dental Waterline Test Kit, HPC Lab Test - 8 vials

**IMS-2008**

Dental Waterline Test Kit, HPC Lab Test - 16 vials

**IMS-2016**

**For more information call 1-800-Hu-Friedy or contact your authorized HuFriedyGroup representative. For question or concerns, contact DUWL Technical Support at [watersupport@hu-friedy.com](mailto:watersupport@hu-friedy.com).**

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted.

©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-516/0423