

Technical Data Sheet

Dental Dams, Non-Latex



Product	Dental Dam, Non-Latex					
Product Overview	Used to isolate and protect teeth and mouth during a dental procedure					
Materials	Manufactured from polyisoprene					
Details	<ul style="list-style-type: none"> • Peppermint flavored • Use of non-latex products drastically reduces risk of allergic reactions for both patients and staff • High tear-strength • 100% Non-chlorinated • Does not contain carbamate esters • Powder-free with less than 5mg/sheet. • Low protein: extractable protein content is below FDA allowable limit of 50µg/sheet (ASTM D5712) 					
Storage Requirements	Keep dry and away from sunlight					
Shelf Life	5 years					
Packaging	<table border="1"> <thead> <tr> <th>Product Code</th> <th>Description</th> </tr> </thead> <tbody> <tr> <td>19500</td> <td>Peppermint flavored, 6 in x 6 in, Medium gauge, 52 Sheets/Box</td> </tr> </tbody> </table>		Product Code	Description	19500	Peppermint flavored, 6 in x 6 in, Medium gauge, 52 Sheets/Box
Product Code	Description					
19500	Peppermint flavored, 6 in x 6 in, Medium gauge, 52 Sheets/Box					
Regulatory Information	Product Code: EIE Classification: Class 1, 510(k) Exempt Regulation: 872.6300, Rubber Dam and Accessories Medical Device Listing: D142786 Made in Malaysia Single-use					