


# SPARKLE™ AIR/WATER SYRINGE TIPS



**Locking indicator and groove** securely installs the tip




**Beveled edge** increases patient comfort



**Quad Core™ Air Channel** ensures air/water separation



**Rounded non-scratch end** prevents o-ring damage



Available in  
6 colors

| SKU             | Description                                       | Quantity |
|-----------------|---|----------|
| <b>BCSAWSCL</b> | Clear   | 250/Bag  |
| <b>BCSAWSWH</b> | White<br>(opaque, not translucent.)               | 250/Bag  |
| <b>BCSAWS</b>   | Assorted<br>(Blue, Purple, Yellow and Amber)      | 250/Bag  |
| <b>BCSAWSB</b>  | Assorted Bulk<br>(Blue, Purple, Yellow and Amber) | 1500/Bag |

## Reduce patient risk and increase workflow productivity.

Studies have shown that the irregular inner surfaces of traditional metal air/water syringe tips allow for the accumulation of deposits and possibly bioburden.<sup>1</sup> This poses a challenge for reliable sterilization as the inside channel is hard to reach and cannot be visually inspected; potentially putting patients at risk for cross-contamination. Single-use Sparkle™ Air/Water Syringe Tips help eliminate the risks associated with reusable syringe tips as well as improve daily workflow and patient comfort.

- **Save time:** No need to reprocess metal syringe tips between patients.
- **Single-use:** Disposable tips eliminate the risk of cross-contamination.
- **No adapter required:** Locking groove ensures tip is securely installed.
- **Efficient:** Ensures distinct air/water separation for every procedure.
- **Durable:** Rigid construction prevents bending during cheek retraction.
- **Translucent colors:** Allow for visualization of water flow and/or blockages.

Learn more at [HuFriedyGroup.com](https://www.HuFriedyGroup.com)

1) Puttaiah, P., Cottone, J., Guildersleeve, J., Azmoudeh, A., Tenney, J. "Rational for Using Single-Use Disposable Air/Water Syringe Tips." Compendium, November 1999, Vol. 20 No.11:1-7.

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | [HuFriedyGroup.com](https://www.HuFriedyGroup.com)

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted.

©2022 Hu-Friedy Mfg. Co., LLC. All rights reserved. T01290 Rev A 0322

Crosstex is now a proud member of

**HuFriedyGroup**  
The Best In Practice