



HU-FRIEDY KEY OPINION LEADER

DR. CORTELLINI AND PROF. TONETTI
MICROSURGERY KIT



Dr. Pierpaolo Cortellini



Prof. Maurizio S. Tonetti

The primary goal of TangramOdis srl is to enhance the knowledge, skills and practice of dental clinicians and dental hygienists.

Our Continuing Education Program includes lectures, group discussions, hands-on exercises with specially developed simulators and unique live 3D surgical demonstrations.

Ample time is allotted for participants to interact with the lecturers and the instructors to allow immediate assimilation of this new information into your practice.

www.tangramodis.it*



“ Composed of essential instruments, each one chosen for a precise purpose, in the treatment of soft tissue in microsurgery. ”

–DR. CORTELLINI AND PROF. TONETTI

MICROSURGERY KIT

CHOOSE THE COLOR.

Just add the color coding number to the last part of the partcode (replace “_” with the number)

- 0 = Red
- 1 = Grey
- 2 = White
- 3 = Orange
- 4 = Purple
- 5 = Yellow
- 6 = Lavander
- 7 = Black
- 8 = Blue
- 9 = Green
- 11 = Ocean Blue

IMEDIN18_S
IMS Infinity
Cassette DIN,
18 Instruments
Surgery



HU-FRIEDY IMS INSTRUMENT MANAGEMENT SYSTEM

By standardising your instruments with IMS, additional time and revenue is generated while extending the useful life of your Hu-Friedy instruments. Dental professionals state that by implementing IMS, they save 5 to 10 minutes** per procedure allowing the staff to spend more time on value-added activities. In addition, safer transportation of soiled instruments reduces risk of injury and exposure to blood borne pathogens. **When compared to single instrument decontamination



PCPUNC156


Probe UNC #UNC15 hdl #6
Qulix1-2-3..14-15



P2N6

Probe Nabers #2N
Bifurcation

**DR. CORTELLINI AND PROF. TONETTI
MICROSURGERY KIT**



SHDPV

Micro Surgical Scalpel Handle
hdl #SinusLine

10-130-05E

Scalpel Handle #5E,
European handle

K012KPO3A6

Knife Gingivectomy Allen #1/2
hdl #6 thin

PPLABANCA6

Periosteal Labanca #3 hdl #6

PPAELA

Periosteal Allen hdl #6 anterior

TPDAPV

Tissue Pliers Micro hdl# SinusLine
straight, diamond powdered, 18 cm

SPV

Scissors Micro Surgical
hdl #SinusLine, 18 cm
ONLY FOR SOFT TISSUE!



10-199-00

Scalpel Blade Remover

MM4

Micro Surgical Mirror round
3mm, flexible

MM7

Micro Surgical Mirror rectangular
7x2 mm, flexible

NHDPV

Needle Holder Micro Surgical
Velvart, straight, diamond dusted,
6-0 - 8-0, 18 cm

S18

Scissors Iris #18 fine
curved 11,5 cm

CURETTES MICRO MINI

POINTS OF PERFORMANCE

- Longer terminal shank designed to access deep periodontal pockets
- Blade is 20 % thinner than a Mini Five Gracey to help reduce tissue distension and ease gingival insertion



IME14DIN5_

IMS Infinity
Cassette 1/4 DIN
for 5 Instruments

Standard
vs. Micro
Mini Five



SMS1/291E2

Curette Gracey
Micro #1/2
hdl #9 MF
oral/labial,
grey

SMS7/897E2

Curette Gracey
Micro #7/8
hdl #9 MF
oral/buccal,
green

SMS11/1293E2

Curette Gracey
Micro #11/12
hdl #9 MF
mesial,
orange

SMS13/1498E2

Curette Gracey
Micro #13/14
hdl #9 MF
distal,
blue

Manufacturer: Hu-Friedy Mfg. Co., LLC | 3232 N. Rockwell Street | Chicago, IL 60618 | USA
Customer Care Department: Hu-Friedy Mfg. Co., LLC.
Free Call +800 HUFRIEDY (00800 48 37 43 39) | Free Fax 00800 48 37 43 40
E-Mail: info@hu-friedy.eu | Website: Hu-Friedy.eu

©2019 Hu-Friedy Mfg. Co., LLC. All rights reserved. HF-GB-CTMS/0819

How the best perform

