

# Surgical Kit

## Luxating Hybrid Kit, Black Line (EPTKITX)



- 1** **IMN5081**  
**8 Instrument Infinity Series™ Cassette, Gray**  
IMS™ Cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



- 2** **EPTLGMX**  
**Luxating Hybrid - Large Mesial, Black Line**  
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- 3** **EPTLGDX**  
**Luxating Hybrid - Large Distal, Black Line**  
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- 4** **EPTSMMX**  
**Luxating Hybrid - Small Mesial, Black Line**  
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- 5** **EPTSMDX**  
**Luxating Hybrid - Small Distal, Black Line**  
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- 6** **EPTSMCX**  
**Luxating Hybrid - Small Curved, Black Line**  
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- 7** **EPTSMSX**  
**Luxating Hybrid Elevator, Small Straight, Black Line**  
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- 8** **EPTSSPX**  
**Luxating Hybrid Elevator, Straight Spade, Black Line**  
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.

**Luxating Hybrid Kit, Black Line (EPTKITX)**



Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | [HuFriedyGroup.com](http://HuFriedyGroup.com)  
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.  
©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-342/1023