

Surgical Kit

Bone Grafting Kit (IMBGKIT)



- 1 IMN4128**
12 Instrument Signature Series™ Cassette, Blue
 IMS™ Cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



- 2 10-130-05**
5 Scalpel Handle, Round Straight
 Used for cutting gingival tissue and making surgical incisions. Provides greater rotational movement with fingertip control and enhances dexterity. Handle measures 6" long. Compatible with standard scalpel blade sizes.



- 3 CSEUWES6**
Euwe Bone Scoop/Plugger, #6 Handle
 Designed for easy transport and non-stick compaction of bone graft material. 3.0mm/7.0mm.



- 4 SP20**
20 Corn Suture Pliers
 Ideal for use with guided tissue regeneration procedures. 15cm



- 5 P96**
9 Molt Periosteal Elevator
 Designed for reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. 4.3mm/8.4mm



- 6 PPBUSER**
Buser Periosteal, #6 Handle
 Designed for reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. 3.8mm/4.3mm



- 7 PPBUSE1/36**
PJO/1 & PJO/3 Combined 13KL, #6 Handle
 Used to remove autologous bone for transplantation material with 4mm and 6mm wide working ends.



- 8 PPBUSE2/36**
PJO/2 & PJO/3 Combined 13KL, #6 Handle
 Used to remove autologous bone for transplantation material with 5mm and 6mm wide working ends.



- 9 CL866**
86 Lucas DE Surgical Curette, #6 Handle
 Designed for curettage, cyst removal and tooth socket debridement. 2.8mm



- 10 SPR1/26**
PRI/2 Prichard Curette, #6 Handle
 Larger and heavier curette for the removal of granulation tissue and tenacious subgingival deposits measuring 7.5" long with 1.5mm wide working ends.



- 11 PMPI**
Membrane Placement Instrument
 Micropoint tip for precise placement and manipulation of membranes.




- 12 MCUPE**
Measuring Cup 25ml with Marking



- 13 S16**
16 Goldman-Fox Scissors
 Used for trimming tissue or cutting sutures. 12.5cm with 1 serrated blade.




- 14 NHHB**
Hegar-Baumgartner Needle Holder
 Used to grasp and guide the needle when suturing. 14 cm. Recommended suture sizes: 1-0, 2-0, 3-0, 4-0



- IM1002**
6-Pack Hinged Instrument Clip
 Secure and maintain open position of scissors, hemostats, and other hinged instruments in lid of cassette.

Bone Grafting Kit (IMBGKIT)



Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.
©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-337 /1023