

# Surgical Kit

## Lateral Sinus Lift Kit (IMLATKIT)



- 1 IM4162**  
**16 Instrument Signature Series™ Cassette, White**  
IMS™ Cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



- 2 IMP65786**  
**Kramer-Nevins Sinus Lift Instrument, Obtuse Angles, #6 Handle**  
Used during a sinus lift procedure to separate/reflect the Schneiderian membrane from the maxillary bone and to elevate the membrane. With obtuse angles 4.1mm/7.6mm.



- 3 IMP6577SC6**  
**Kramer-Nevins Sinus Lift Instrument, Acute Angles, #6 Handle**  
Used during a sinus lift procedure to separate/reflect the Schneiderian membrane from the maxillary bone and to elevate the membrane. With acute angles 4.1mm/7.6mm.



- 4 BPACK**  
**Bone Packer DE 6-8mm**  
Used to condense bone during grafting procedures with working ends 6mm and 8mm wide.



- 5 SYBG**  
**Syringe Bone Grafting 7mm**  
Used during grafting procedures. Inner Diameter 5.8 mm. Load from tip to prevent clogging.



- 6 BONT**  
**Bone Tamper 7mm**  
Used to separate the Schneiderian membrane from surrounding bone during lateral sinus lift procedures.



- 7 NHCW**  
**Crile-Wood Needle Holder**  
Used to separate the Schneiderian membrane from surrounding bone during lateral sinus lift procedures.



- 8 S16**  
**16 Goldman-Fox Scissors**  
Used for trimming tissue or cutting sutures. 12.5cm with 1 serrated blade.



- 9 CRURBAN**  
**Urban Cheek Retractor**  
Used to hold mucoperiosteal flaps, cheeks, lips and tongue away from the surgical area. 20.6mm



- 10 IM1002**  
**6-Pack Hinged Instrument Clip**  
Secure and maintain open position of scissors, hemostats, and other hinged instruments in lid of cassette.

**Lateral Sinus Lift Kit (IMLATKIT)**



Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com  
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.  
©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-340 /1023