

# Surgical Kit

## Microsurgical Kit, Black Line (IMNMICSURX)



- 1 IMN4121**  
**12 Instrument Infinity Series™ Cassette, Gray**  
IMS™ cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



- 2 10-130-70**  
**7 Mini Scalpel Blade Holder**  
Designed to hold micro scalpel blades for minimally invasive procedures and microsurgery.



- 3 S13K/13KL6**  
**13K/13KL DE Surgical Chisel**  
Used to recontour periodontal bone or to remove remaining tooth fibers during periodontal surgery with curved working ends 3.2mm and 4.9mm wide.



- 4 C36/37**  
**36/37 Back-Action DE Chisel**  
Designed for removing bone adjacent to the tooth without causing trauma and is especially useful on the distal of last molars with working ends 3.8mm wide.



- 5 8-905DD**  
**Precision Dressing Forceps, Diamond Dusted**  
Used to grasp and stabilize soft tissue, 7" length.



- 6 NH5024CORT**  
**Cortellini Micro Castro Needle Holder**  
Used to grasp and guide the needle when suturing.



- 7 PPBUSERX**  
**Buser Periosteal, Satin #6 Handle**  
Featuring a sharp round blade and narrow pointed end for interproximal use. For reflecting and retracting the mucoperiosteum.



- 8 SPVX**  
**Microsurgical Scissors**  
Used for trimming tissue or cutting sutures. 16.5cm.



- 9 PCPUNC156**  
**CC SE UNC Probe, 1-15, #6 Handle**  
Single-ended, color-coded probe. ergonomic, Satin Steel™ handle. University of North Carolina.



- 10 MIR5**  
**5 FS CS Mouth Mirror**  
Single-sided Cone Socket mirror. 15/16 in. (24mm) diameter.



- 11 MH6**  
**Cone Socket Mirror Handle, Satin Steel**  
Used with Single and Double-Sided Cone Socket Mirrors. The blunt end of the mirror handle is ideal for testing sensitivity to percussion and assessing tooth mobility.

Microsurgical Kit, Black Line (IMNMICSURX)



Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com  
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.  
©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-349/1023