

Surgical Kit

RESCUE Kit (IMRESCUE)



- 1 IM4160**
16 Instrument Signature Series™
Cassette, Red
IMS™ Cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



- 2 IMS-1372SH**
IMS 1/2 Bur Cushion - Short Lid
Designed to hold burs securely in place during cleaning and sterilization. Hollow design allows drainage and proper cleaning.



- 3 EPTSSP**
Luxating Hybrid Elevator, Straight Spade
A unique cross between a periostome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- 4 EPTSMC**
Luxating Hybrid, Small Curved
A unique cross between a periostome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- 5 EPTSMS**
Luxating Hybrid Elevator, Small Straight
A unique cross between a periostome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- 6 CL866**
86 Lucas DE Surgical Curette #6 Handle
Designed for curettage, cyst removal and tooth socket debridement. 2.8mm.




- 7 EHB13/146**
13/14 Heidbrink DE Root Pick
Used for the removal of root tips in difficult-to-reach areas. Designed for delicate root teasing procedures.



- 8 RMF**
Mini Friedman Rongeurs
Used for trimming and recontouring alveolar bone and gross tissue removal.

Rescue Kit (IMRESCUE)



Not pictured: #2 [IMS-1372SH](#)

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | [HuFriedyGroup.com](#)
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.
©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-351/1023