


STAINLESS STEEL MOLAR CROWNS



Experience the perfect fit for all patients. No matter what patient walks through your door, HuFriedyGroup now has you covered more than ever before. **Our line of Stainless Steel Crowns has expanded to include Permanent Molar Crowns!** All of our Stainless Steel Crowns deliver all the features and benefits of our Pedo Crowns that you love, but now in the form of Permanent Crowns as well.



SKU	Description	Quantity
SSC-KIT	Primary Stainless Steel Crowns Introductory Kit	48 Crowns (1 of each size)
SSC-PKIT	Permanent Stainless Steel Crowns Introductory Kit	24 Crowns (1 of each size)
SSC-PKITM	Permanent Stainless Steel Crowns Master Kit	96 Crowns (4 of each size)

Pre-trimmed, crimped and contoured surfaces for less time spent chairside

Accurate occlusal anatomy mimicking the natural tooth

Ideal wall height meaning less trimming and crimping (if at all!)




Optimal occlusal thickness with superior resilience to abrasion and perforation


More narrow mesial-distal maxillary first primary molar width for quicker adjacent crown placement

Stainless Steel Crowns


Refill Part Codes




SSC-ULD7 SSC-ULD6 SSC-ULD5 SSC-ULD4 SSC-ULD3 SSC-ULD2



SSC-URD2 SSC-URD3 SSC-URD4 SSC-URD5 SSC-URD6 SSC-URD7



SSC-ULE7 SSC-ULE6 SSC-ULE5 SSC-ULE4 SSC-ULE3 SSC-ULE2





SSC-URE2 SSC-URE3 SSC-URE4 SSC-URE5 SSC-URE6 SSC-URE7




SSC-6UL7 SSC-6UL6 SSC-6UL5 SSC-6UL4 SSC-6UL3 SSC-6UL2




SSC-6UR2 SSC-6UR3 SSC-6UR4 SSC-6UR5 SSC-6UR6 SSC-6UR7




SSC-LLD7 SSC-LLD6 SSC-LLD5 SSC-LLD4 SSC-LLD3 SSC-LLD2



SSC-LRD2 SSC-LRD3 SSC-LRD4 SSC-LRD5 SSC-LRD6 SSC-LRD7



SSC-LLE7 SSC-LLE6 SSC-LLE5 SSC-LLE4 SSC-LLE3 SSC-LLE2



SSC-LRE2 SSC-LRE3 SSC-LRE4 SSC-LRE5 SSC-LRE6 SSC-LRE7



SSC-6LL7 SSC-6LL6 SSC-6LL5 SSC-6LL4 SSC-6LL3 SSC-6LL2



SSC-6LR2 SSC-6LR3 SSC-6LR4 SSC-6LR5 SSC-6LR6 SSC-6LR7



For more information, visit [HuFriedyGroup.com/Crowns](https://www.hufriedygroup.com/Crowns)