

# PERMA SHARP® NEEDLE HOLDERS

by Hu-Friedy

## FACT


Clinicians require needle holders that are reliable and provide a firm, positive needle grip to suture securely.

## CHALLENGE

Needle holders that do not provide optimal needle engagement increase the chance of rotation and or slippage.

## SOLUTION

Hu-Friedy's Perma Sharp needle holders feature tungsten carbide inserts with a matched diamond-cut pyramidal profile for enhanced longevity and a secure, positive grip.



## POINTS OF PERFORMANCE

- **Tungsten carbide inserts** with a diamond cut pyramidal profile ensure a firm positive needle grip
- **Reduced risk of cutting the suture material** with specially bonded and matched carbide inserts
- **Optimal grip and control** with a variety of design options, including standard, finger and palm controlled closure
- **Longevity and corrosion resistance** afforded by uncompromising quality



How the best perform

**Hu-Friedy**

## NEEDLE HOLDER REFERENCE CHART

HU-FRIEDY PART CODE	DESCRIPTION	NEEDLE HOLDER SIZE	RECOMMENDED SUTURE SIZES
NH5022	Arruga	6-1/4" (16 cm)	5-0, 6-0, 7-0, 8-0
NH5030	Webster Smooth*	4-1/2" (11.5 cm)	8-0, 9-0, 10-0
NH5032	Derf	4-1/2" (11.5 cm)	3-0, 4-0, 5-0
NH5034	Hegar-Baumgartner	5-1/2" (14 cm)	1-0, 2-0, 3-0, 4-0
NH5036	Halsey Smooth*	5" (13 cm)	8-0, 9-0, 10-0
NH5037	Halsey	5" (13 cm)	3-0, 4-0, 5-0, 6-0
NH5038	Crile-Wood	6" (15 cm)	3-0, 4-0, 5-0, 6-0
NH5042	Mayo-Hegar	6-1/4" (16 cm)	1-0, 2-0, 3-0, 4-0
NH5044	Mayo-Hegar	7" (18 cm)	1-0, 2-0, 3-0, 4-0
NH5046	DeBakey Delicate	7" (18 cm)	4-0, 5-0, 6-0
NH5050	HMT Thin Jaw	5-1/2" (14 cm)	4-0, 5-0, 6-0
NH5052	Fine Swedish Pattern	6" (15 cm)	5-0, 6-0, 7-0, 8-0
NH5068	Olsen-Hegar Needle Holder Scissors	6-3/4" (17 cm)	3-0, 4-0, 5-0, 6-0
NH5072	Olsen-Hegar Needle Holder Scissors	5-1/2" (14 cm)	3-0, 4-0, 5-0, 6-0
NH5074	Mathieu-Kocher	5-1/2" (14 cm)	1-0, 2-0, 3-0, 4-0
NH5076	Mathieu	5-1/2" (14 cm)	1-0, 2-0, 3-0, 4-0
NH5082	Micro-Vascular Pattern	6" (15 cm)	6-0, 7-0, 8-0
NH5084	Micro-Vascular Pattern	7" (18 cm)	6-0, 7-0, 8-0
NH5092	Mini-Ryder Pattern	6" (15 cm)	8-0, 9-0, 10-0
NH5094	Mini-Ryder Pattern	7" (18 cm)	8-0, 9-0, 10-0
NH5020	Castroviejo Straight	5-1/2" (14 cm)	5-0, 6-0, 7-0
NH5020R	Castroviejo Straight, Round Handle	5-1/2" (14 cm)	5-0, 6-0, 7-0
NH5021	Castroviejo Curved	5-1/2" (14 cm)	5-0, 6-0, 7-0
NH5024	Castroviejo Straight	7" (18 cm)	5-0, 6-0, 7-0, 8-0, 9-0, 10-0
NH5024R	Castroviejo Straight, Round Handle	7" (18 cm)	5-0, 6-0, 7-0, 8-0, 9-0, 10-0
NH5024RC	Castroviejo Curved, Round Handle	7" (18 cm)	5-0, 6-0, 7-0, 8-0, 9-0, 10-0
NH5024SIM	Simion Castroviejo Straight, Round Handle	7" (18 cm)	5-0, 6-0, 7-0, 8-0, 9-0, 10-0
NH5020M	Micro Castroviejo Straight	5-1/2" (14 cm)	5-0, 6-0, 7-0
NH5021M	Micro Castroviejo Curved	5-1/2" (14 cm)	5-0, 6-0, 7-0
NH5024CORT	Cortellini Micro Castroviejo Straight	7" (18 cm)	5-0, 6-0, 7-0, 8-0
NHM-5025R	Micro Castroviejo Straight, Diamond Dusted	6" (15 cm)	5-0, 6-0, 7-0, 8-0, 9-0, 10-0
NHM-5026R	Micro Castroviejo Straight, Diamond Dusted	7" (18 cm)	6-0, 7-0, 8-0, 9-0

\* Non-Serrated

Hu-Friedy also offers an array of standard needle holders.

Ask your authorized Hu-Friedy dealer for details, or call us at **1-800-HU-FRIEDY**.

Visit us online at **Hu-Friedy.com**

Call 1-800-Hu-Friedy or contact your authorized Hu-Friedy representative for more information.  
©2013 Hu-Friedy Mfg. Co., LLC. All rights reserved. HF-250/0513

