


PART CODE	DESCRIPTION
IMP204SDT	204SD
IMPLG1/2T	Langer 1/2
IMPM1/2T	Mini Five 1/2
IMPM11/12T	Mini Five 11/12
IMPM13/14T	Mini Five 13/14

TITANIUM IMPLANT SCALERS

PERFORM AT YOUR BEST with Hu-Friedy's new, state of the art line of Titanium Implant Scalers. Made from the same titanium alloy as implants and abutments, these scalers are expertly designed and manufactured for implant maintenance, debridement, biofilm removal and can be used both supra and subgingival.

- **Improved visual acuity** with Hu-Friedy's signature Teal colored Anodized Titanium
- **Optimized sharpness** with cutting edges that are finely honed and sharpened post Anodization process
- **Enhanced instrument characteristics** for Implant maintenance with the ability to remove calculus and biofilm
- **Increased instrument value** with the ability to be sharpened at any time to extend the useful life of the instrument
- **Ergonomic sensitivity** featuring our Ergo handle that offers a large diameter and cross/hatch design for maximum tactile grip and support



VISIT US ONLINE AT HU-FRIEDY.COM/TIS

©2017 Hu-Friedy Mfg. Co., LLC. All rights reserved. [820] N2670/0817

How the best perform

