

PEDIATRIC FORCEPS



PRECISION DESIGNS FOR
CHILD-FRIENDLY EXTRACTIONS



Hu-Friedy



Hu-Friedy's complete line of pediatric forceps provides child-friendly extraction solutions to enhance patient comfort and preserve delicate tissue and bone. These are precision-designed specifically for a child's anatomical tooth structure, ensuring clinical performance and optimal pediatric dental care.

Crafted from the finest materials, Hu-Friedy's extensive line of pedo forceps include several options to suit your clinical preference.

LIFETIME WARRANTY

Call 1-800-Hu-Friedy or visit us online at Hu-Friedy.com for more information about the Presidential and Rainbow Presidential Pedo Series Lifetime Warranty.*

STANDARD PEDO SERIES

The precise, anatomical design assures a firm, accurate grasp of primary tooth and root surfaces for smooth, seamless extractions with a reduced risk of complications or slipping. Hand-crafted from surgical-grade stainless steel, Hu-Friedy's pediatric forceps resist corrosion to provide superior performance and longevity.

PRESIDENTIAL® PEDO SERIES

Similar to the Standard Pediatric Series, the Hu-Friedy Presidential Forceps are meticulously crafted and finished with a precision hinge, carrying a lifetime warranty against defects in material and workmanship.


RAINBOW PRESIDENTIAL PEDO SERIES

Hu-Friedy's Rainbow Presidential Pediatric Forceps feature a unique rainbow of vibrant colors to reduce child anxiety creating a more positive patient experience. The series carries a lifetime warranty against material and workmanship defects.*

KINDER PEDO SERIES

The Hu-Friedy Kinder Pediatric Forceps fit easily and comfortably in your hand. Smaller than the Standard Pediatric Forceps, these are easily hidden from a child's view. Created for small, shallow crowns and conical shaped teeth, the Kinder design requires less force, decreasing the risk of tissue damage and broken crowns.

*The special rainbow finish on the F150SR and F151SR must be handled with care to avoid scratching and chipping. These forceps will not be repaired or replaced due to scratches or chips to the finish unless there is a manufacturing defect present or the product has not been used.





ORDERING INFORMATION

STANDARD PEDO SERIES

F1S	#1 Pedo Forceps
F16S	#16S Pedo Forceps
F23S	#23S Pedo Cowhorn Forceps
F39	#39 Pedo Forceps
F40	#40 Pedo Forceps
F150XAS	#150S Pedo Forceps Upper Split Beaks Serrated
F151XAS	#151S Pedo Forceps Lower Split Beaks Serrated

PRESIDENTIAL PEDO SERIES

F150S	#150S Pedo Forceps Upper
F151S	#151S Pedo Forceps Lower

RAINBOW PRESIDENTIAL PEDO SERIES

F150SR	#150SR Pedo Forceps Rainbow Upper
F151SR	#151SR Pedo Forceps Rainbow Lower

KINDER PEDO SERIES

F23K	#23K Pedo Forceps Cowhorn
F150K	#150K Pedo Forceps Upper
F151K	#151K Pedo Forceps Lower



Hu-Friedy Mfg. Co., LLC
3232 N. Rockwell St.
Chicago, IL 60618

1-800-Hu-Friedy
Hu-Friedy.com