



Contains One

Air/Water Syringe Clip

for IMS® Resin Cassettes


Part Code: IMS-1278

Stainless Steel

An **IMS**® Product **CLEAN AND STERILIZE BEFORE USE.**

CAUTION: Overextending, bending or attempting to modify can cause breakage. Clip is recommended for use with large IMS Resin Cassettes: IMS-1112, IMS-1116, IMS-1118, IMS-1124. Also for use with small open ended IMS Resin Cassettes: IMS-3107, IMS-3112, IMS-3116.

INSTRUCTIONS



Place the flat side on the outside of the cassette.



MADE IN USA

Manufactured for:
Hu-Friedy Mfg. Co., LLC
3232 N. Rockwell St., Chicago, IL 60618 USA
1-800-Hu-Friedy | Hu-Friedy.com

Hu-Friedy Mfg. Co., LLC
Zweigniederlassung Deutschland
Ziegeleiweg 1
D-78532 Tuttlingen-Möhringen
A4435/1010