

INSTRUMENT  
MANAGEMENT

**AIR/WATER SYRINGE  
CASSETTE CLIP**

Part Code: IM1005 and IM1005N



CONTAINS ONE  
CLEAN AND STERILIZE BEFORE USE.

***Hu-Friedy***<sup>®</sup>

# INSTRUCTIONS

Use clip on side of cassette. Insert air/water syringe tip into clip.

**CAUTION: Overextending, bending or attempting to modify can cause breakage. Clip is recommended for use with IMS Signature Series Cassettes only.**

**Manufactured for Hu-Friedy Mfg. Co., LLC**  
3232 N. Rockwell St. | Chicago, IL 60618 | USA  
1-800-Hu-Friedy | Hu-Friedy.com

**Hu-Friedy Mfg. Co., LLC**  
Zweigniederlassung Deutschland  
Ziegeleiweg 1  
D-78532 Tuttlingen-Möhringen

