



Contains One

### Air/Water Syringe Clip

for IMS® Resin Cassettes


Part Code: IMS-1279

Stainless Steel

An **IMS**® Product **CLEAN AND STERILIZE BEFORE USE.**

CAUTION: Overextending, bending or attempting to modify can cause breakage. Clip is longer than IMS-1278 for proper fit on small IMS Resin Cassettes with closed sides: IMS-3107, IMS-3112, IMS-3116.

### INSTRUCTIONS



Place the flat side on the outer (closed side) of the resin cassette.



**MADE IN USA**

Manufactured for:  
Hu-Friedy Mfg. Co., LLC  
3232 N. Rockwell St., Chicago, IL 60618 USA  
1-800-Hu-Friedy | Hu-Friedy.com

Hu-Friedy Mfg. Co., LLC  
Zweigniederlassung Deutschland  
Ziegeleiweg 1  
D-78532 Tuttlingen-Möhringen  
B7847/0211