

Slim Flush Cut & Hold Distal End Cutter
GET **CLOSER** TO A BETTER RESULT



Hu-Friedy[®]

IMPROVED ACCESS. GREATER VISIBILITY.

With the Slim Flush Cut & Hold instrument, Hu-Friedy has combined artistry and function to produce a product that is uniquely capable of adding value to the Orthodontist's ability to deliver optimal patient care. Some of the key benefits include:

- **Improved access and enhanced visibility** made possible by a working end that has been reduced by 50% as compared to 678-111
- **Optimized patient comfort and safety and significant time savings,** delivered through a true flush cut & hold instrument that helps decrease emergency stick appointments and helps improve efficiency by eliminating the need for a cotton roll or piece of gauze to catch the cut portion
- **Additional flexibility on archwire selection** given the instrument's ability to cut up to .021" x .025"
- **Steam sterilizable**

SLIMMER BY COMPARISON – 50% THINNER END


678-111
Flush Cut & Hold
Distal End Cutter




678-113
Slim Flush Cut & Hold
Distal End Cutter



REACH FOR A REMARKABLY FLUSH CUT.



Buccal view of the Slim Flush Cut & Hold Distal End Cutter cutting the archwire flush to the buccal tube.



Occlusal view of the Slim Flush Cut & Hold Distal End Cutter cutting the archwire flush to the buccal tube.

INSTRUMENT PART-CODES



678-113
Slim Flush Cut & Hold
Distal End Cutter



678-111
Flush Cut & Hold
Distal End Cutter

Accessories Part-codes:

| | |
|---------|---|
| 678-013 | Silicone Insert for 678-113 (pack of 3) |
| 678-011 | Silicone Insert for 678-111 (pack of 3) |

FREE SHARPENING OFFER

All Hu-Friedy cutters qualify
for one free sharpening within
12 months of purchase*

Instrument Check-Up Guide

Check all your instruments for the following:

- ☐ Dull cutting edge
- ☐ Broken working end
- ☐ Bent or misshaped working end
- ☐ Corrosion
- ☐ Pitting

To view product in action visit Hu-Friedy.com. Call 1-800-Hu-Friedy
or contact your Hu-Friedy representative for more information.



Hu-Friedy Mfg. Co., Inc.
3232 N. Rockwell St.
Chicago, IL 60618
1-800-Hu-Friedy
Hu-Friedy.com