

## IMS HINGED INSTRUMENT CLIP

Part Code: IM1000 | Stainless Steel



CONTAINS TWO  
CLEAN AND STERILIZE BEFORE USE.



### INSTRUCTIONS


Use clip on lid of cassette.

**CAUTION:** Overextending, bending or attempting to modify can cause breakage. Clip is recommended for use with IMS Infinity & Signature Series Cassettes only.

**Manufactured for Hu-Friedy Mfg. Co., LLC**  
3232 N. Rockwell St. | Chicago, IL 60618 | USA  
1-800-Hu-Friedy | Hu-Friedy.com

**Hu-Friedy Mfg. Co., LLC**  
Zweigniederlassung Deutschland  
Kleines Öschle 8, D-78532 Tuttlingen

MADE IN USA



B2932/0617