


Disposable Prophy Angles

Twist™ Plus Prophy Angles



POINTS OF PERFORMANCE

- The Twist™ Plus Prophy Angle is the ONLY disposable that oscillates 120° in a smooth, continuous counter-rotational motion to eliminate the splatter and frictional heat associated with traditional 360° spinning prophy angles.
- Single-use design for infection control compliance.
- Round head design improves visibility inside the oral cavity and increases maneuverability and patient comfort.
- Superior cup design with external fluting: Enhanced inner web provides better pickup and retention of paste, while external fluting improves interproximal stain/biofilm removal.
- Recessed ridges in the grip makes a modified pen grasp more comfortable and secure while wearing gloves.
- Handpiece connection is on the underside, which reduces the potential for hand movements to push the angle off the handpiece.
- Tight seal of the cup to head reduces the risk of catching soft tissue in the corners of the mouth.

SKU	Description	Quantity
TPLUSPASC	Twist™ Plus DPA, Soft Cup	100/box
TPLUSPAFC	Twist™ Plus DPA, Firm Cup	100/box
TPLUSPATB	Twist™ Plus DPA, Tapered Brush	100/box

Disposable Propphy Angles

Twist™ Propphy Angles




POINTS OF PERFORMANCE

- Ergonomically designed angle that provides a 90° reciprocating “back and forth” rotating motion for continuous contact.
- Smooth back and forth motion reduces frictional heat, keeping the tooth surface cooler and reduces splatter.
- Soft web cup and ribbed flair edges, provides easy application and spreading of prophyl paste.
- Single-use.
- Propphy angles not made with natural rubber latex.
- Made outside of the USA.

SKU	Description	Quantity
TPAFC	Twist™ Propphy Angle, Firm	100/box
TPASC	Twist™ Propphy Angle, Soft	100/box
TPATB	Twist™ Propphy Angle, Tapered Brush	100/box
TPAFB	Twist™ Propphy Angle, Flat Brush	100/box
TPASCB	Twist™ Propphy Angle, Soft	500/box

Sparkle™ Propphy Angles



POINTS OF PERFORMANCE

Created with hygienists' feedback:

- Patented beveled gear design runs vibration-free/ noise-free.
- Flared cup design improves subgingival access and reduces spatter, minimizing cross-contamination.
- Light colored angle with dark colored cup provides greater visibility.
- Single-use.

SKU	Description	Quantity
TPAFEZ	Sparkle™ EZ Contra Angle, Firm	100/box
TPASEZ	Sparkle™ EZ Contra Angle, Soft	100/box
TPASFC1	Sparkle™ Propphy Angle, Firm	100/box
TPASSC1	Sparkle™ Propphy Angle, Soft	100/box
TPATBC1	Sparkle™ Propphy Angle, Tapered Brush	100/box
TPASFC5	Sparkle™ Propphy Angle, Firm	500/box
TPASSC5	Sparkle™ Propphy Angle, Soft	500/box

For more information, visit [HuFriedyGroup.com](https://www.HuFriedyGroup.com)

Crosstex is now a proud member of