

# SLIM ELEVATORS

Designed to be  
**13-20% thinner**  
than standard  
luxating  
elevators!

**Slim luxating elevators are ideal for atraumatic extractions.**

The slender working end glides along the alveolar bone and severs the periodontal ligaments during surgical extractions.



**ELSLIM2CLX**

**Slim Elevators, Long**

**ELSLIM2CSX**

**Slim Elevators, Short**

## Features & Benefits

- Available in **multiple sizes, lengths, and shapes** providing customized treatment based on patient need
- **Lightweight, large diameter handle** for maximum comfort, reduced hand fatigue and increased control\*
- Available in **Black Line and Standard Stainless-Steel**
  - Performance engineered coating for a harder, smoother surface for **enhanced lubricity and optimal edge retention**
  - Distinct black finish for **enriched contrast** between the instrument, surgical site and underlying tissue
  - **Reduced light reflection** afforded by a matte finished handle and black working ends

Part Code	Description
<b>ELSLIM2CL</b>	Slim Luxating Elevator 2MM, Long Curved
<b>ELSLIM2CLX</b>	Slim Luxating Elevator 2MM, Long Curved, Black Line
<b>ELSLIM2CS</b>	Slim Luxating Elevator 2MM, Short Curved
<b>ELSLIM2CSX</b>	Slim Luxating Elevator 2MM, Short Curved, Black Line
<b>ELSLIM2SL</b>	Slim Luxating Elevator 2MM, Long Straight
<b>ELSLIM2SLX</b>	Slim Luxating Elevator 2MM, Long Straight, Black Line
<b>ELSLIM2SS</b>	Slim Luxating Elevator 2MM, Short Straight
<b>ELSLIM2SSX</b>	Slim Luxating Elevator 2MM, Short Straight, Black Line
<b>ELSLIM3CL</b>	Slim Luxating Elevator 3MM, Long Curved
<b>ELSLIM3CLX</b>	Slim Luxating Elevator 3MM, Long Curved, Black Line
<b>ELSLIM3CS</b>	Slim Luxating Elevator 3MM, Short Curved
<b>ELSLIM3CSX</b>	Slim Luxating Elevator 3MM, Short Curved, Black Line
<b>ELSLIM3SL</b>	Slim Luxating Elevator 3MM, Long Straight
<b>ELSLIM3SLX</b>	Slim Luxating Elevator 3MM, Long Straight, Black Line
<b>ELSLIM3SS</b>	Slim Luxating Elevator 3MM, Short Straight
<b>ELSLIM3SSX</b>	Slim Luxating Elevator 3MM, Short Straight, Black Line

\*Source: Gerwadowski, L.J., McFall, D.B., Stach, D.: Carpal Tunnel Syndrome; Risk Factors and Preventive Strategies for the Dental Hygienist. Journal of Dental Hygiene, February 1992.

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | Hu-Friedy.com

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted.

©2021 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-577/0721