


# Diamond Dusted Forceps

Our **Diamond Dusted Forceps** feature a fine diamond coating which increases resistance between the forceps and tooth surface providing enhanced engagement. Diamond Dusted Forceps are available in an array of patterns which allows for full-mouth coverage.

**A light diamond coating** bonded to longitudinal beak serrations increases grip by 20% compared to traditional forceps.\*



**Thin, tapered beak designs** allow for improved access compared to traditional forceps beaks.



**Matte finish** for additional contrast and reduced light reflection.

**Unique hole design** reduces the weight of the forceps, providing increased comfort and control.

# Diamond Dusted Forceps

are available in the following Atraumair patterns:



## Upper

**F1XSD**     **1 Standard Forceps, Diamond Dusted**  
Upper Incisors

**FAF1XSD**     **1 Apical Forceps, Diamond Dusted**  
Upper Incisors

**F35XSD**     **35 Forceps, Diamond Dusted**  
Upper Canines and Premolars

**FAF35XSD**     **35 Apical Forceps, Diamond Dusted**  
Upper Canines and Premolars

**FX44XSD**     **44 Root Forceps, Diamond Dusted**  
Upper Roots

**FX51XSD**     **51 Root Forceps, Diamond Dusted**  
Upper Roots

**F67XSD**     **67 Apical Forceps, Diamond Dusted**  
Upper Molars

## Lower

**F13XSD**     **13 Forceps, Diamond Dusted**  
Lower Premolars

**F21XSD**     **21 Forceps, Diamond Dusted**  
1st and 2nd Molars

**F36XSD**     **36 Forceps, Diamond Dusted**  
Lower Incisors

**FAF36XSD**     **36 Apical Forceps, Diamond Dusted**  
Lower Premolars and Incisors

**FX45XSD**     **45 Root Forceps, Diamond Dusted**  
Lower Roots, Serrated

**FX73XSD**     **73 Forceps, Diamond Dusted**  
Lower Molars

**FX79XSD**     **79 Forceps, Diamond Dusted**  
Lower Molars

Find your optimal extraction forceps at [HuFriedyGroup.com/Surgical-Forceps](https://www.HuFriedyGroup.com/Surgical-Forceps)

\*Data on file compared to standard Hu-Friedy Forceps.  
Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com  
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted.  
©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-666/0123