
SURGICAL


Introducing Hu-Friedy's New

ATRAUMAIR FORCEPS



Hu-Friedy's new line of Atraumair Forceps are designed for optimal access and enhanced crown and root engagement.

Each forceps is meticulously crafted to the highest quality standard, ensuring superior performance and longevity. Atraumair Forceps are available in 13 unique patterns that provide a complete line for a variety of extraction cases.



ATRAUMAIR POINTS OF PERFORMANCE:

- Enhanced beak geometries and longitudinal beak serrations provide a **firm, secure grasp of the crown and root reducing the risk of crushing or fracturing the tooth**
- Laser marked hinge for **greater visibility when bagged for sterilization**
- Matte finish for **enhanced contrast and reduced light reflection**
- Unique hole design reduces weight of the forceps, providing **increased comfort and control**
- Meticulously crafted from the finest surgical grade stainless steel

How the best perform Hu-Friedy®

ATRAUMAIR FORCEPS



67 Forceps
Upper Molars
| F67XS



1 Standard Forceps
Upper Incisors
| F1XS



1 Apical Forceps
Upper Incisors
| FAF1XS



35 Forceps, Upper
Canines & Premolars
| F35XS

35 Apical Forceps, Upper
Canines & Premolars
| FAF35XS



51 Root Forceps
Upper Roots, Serrated
| FX51XS



36 Forceps
Lower Incisors
& Canines
| F36XS



36 Apical Forceps
Lower Premolars,
Canines & Incisors
| FAF36XS



13 Forceps
Lower Premolars
| F13XS



79 Forceps
Lower Molars
| FX79XS



73 Forceps
Lower Molars
| FX73XS



23 Cowhorn Forceps
Lower Molars
| F23XS



45 Root Forceps
Lower Roots, Serrated
| FX45XS

Visit us online at **Hu-Friedy.com**

©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. HF-248/0115

Hu-Friedy®