# PINK MATTERS® X HUFRIEDYGROUP SURGICAL KIT


## Minimally Invasive Periodontal and Implant Surgery Solution

Pink Matters® was initiated in Europe by a group of clinicians in order to improve both the patient and practitioner's quality of life. "The Pink Matters Group[1] believes in the synergy of skills. Because of this, our selective range of products are designed by practitioners and developed in collaboration with the most recognized and innovative manufacturers and researchers. The Pink Matters core value is based on ethics and high quality standards through science, practice and knowledge.

**Discover more on Pink Matters Group website**
**www.pink-matters.com**



> Minimally invasive periodontal and implant surgery has become the new standard to improve healing processes and patient quality of life. With the development of many innovative surgical procedures the need for specifically designed instruments has become crucial. Pink Matters®, a group of dedicated clinicians, partnered with HuFriedyGroup, the world leader of high quality surgical instruments, to develop a universal surgical kit, including five new instruments. It allows clinicians to perform all modern periodontal and implant surgery procedures.
>
> **Pink Matters Group**



**PMNANOSS**

**PMNANOAS**

**PMEDGES**

**PMSPLITSX**

**PMPRODIGYS**

## VERSATILITY AND EFFICIENCY: PINK MATTERS KEYS FOR SUCCESS

The Pink Matters® X HuFriedyGroup Surgical Kit is the result of the partnership between the Pink Matters group and HuFriedyGroup. The kit includes five highly versatile instruments designed for minimally invasive periodontal and implant surgery.

In the kit, there are five instruments: two periodontal knives (which can also be used as micro-elevators), a chisel, and an innovative periodontal curette. Instruments benefit from EverEdge™ 2.0 Technology which utilizes optimized heat control and enhanced finishing processes in order to create sharp edges needed for effective procedures and continued use. The periosteal is offered with Black Line coating. This dark, hard coating provides enriched contrast at the surgical site and enhanced lubricity.


HuFriedyGroup
The Best In Practice

# THE PINK MATTERS® INSTRUMENTS

WATCH OUR **TUTORIAL VIDEOS** FOR THE PINK MATTERS® INSTRUMENTS








## PMNANOSS

### PM Nanos S (Straight), EverEdge™ 2.0 Perio Knife, Short

Used for supra periosteal dissection apically and laterally during tunneling and soft tissue surgical procedures. This instrument can also be used to elevate thin papilla during partial or full thickness flap surgery and is convenient for papillary preservation technique.








## PMNANOAS

### PM Nanos A (Angulated), EverEdge™ 2.0 Perio Knife, Short

Used for supra periosteal dissection apically and laterally during tunneling and soft tissue surgical procedures. This instrument can also help in the creation of a partial or full thickness envelope flap during connective tissue grafting procedures in perio plastic surgery and implant surgery. Nano A is also useful when elevating a full thickness flap in the lingual and palatal area during flap surgery.





## PMEDGES

### PM Edge, EverEdge™ 2.0 Crestal Bone Chisel, Short

Designed to elevate full thickness flap in the posterior area in preparation for implant placement. The Edge can be also very useful early on in the implant placement procedure when reflecting the thicker buccal soft tissue. This instrument allows for bone recontouring during crown lengthening and is very efficient at bone harvesting in areas with delicate anatomy, such as the nasal spine.







## PMSPLITSX

### PM Split, Black Line Periosteal, Short

Two sides thoughtfully developed: one side is sharp to safely elevate full thickness flaps in perio and implant surgery. The other side is blunt to reduce the surgical risk when separating elastic fibers during buccal and lingual flap advancement.





## PMPRODIGYS

### PM Prodigy, EverEdge™ 2.0 Surgical Curette, Short

Developed to improve granulation tissue removal during flap surgery. The angulation of the curette improves the access to narrow periodontal and peri-implant defects, particularly in the posterior area.

# ADDITIONAL KIT INSTRUMENTS



**S204SCCHE2**

S204S
Matte finish



**PCPUNC15CH**

Probe UNC
#UNC15 CH
Hdl Qulix, 1-15



**10-130-5EM**

European Round
Scalpel Handle
matte finish



**MIR5DSHDB**

Double-sided Mirror
**Black Line**

**MH6B**

Mirror Handle
**Black Line**

## REDUCES GLARE UP TO 80%[2]

See your patients – not your instruments. Glare interferes with visibility. Reducing this disruption to your visibility allows you to see more detail and offer more accurate patient care. HD Black Line Mirrors are engineered to optimize clinical outcomes by delivering superior visibility throughout any dental procedure.

**SPV**

Curved Microsurgical
Scissor



**S18SC**

Iris Curved Super-Cut™
Scissor



**CRMX**

Minnesota Retractor,
**Black Line**




## FEATURES AND BENEFITS OF BLACK LINE:

- Performance engineered coating for a harder, smoother surface for optimal edge retention and enhanced lubricity
- Distinct black finish for enriched contrast and visual acuity at the surgical site and underlying tissue
- Reduced light reflection afforded by a matte finished handle and black working ends
- Unique, smooth, large diameter lightweight handle for maximum comfort, reduced hand fatigue and increased control
- Handcrafted with Immunity Steel™ Stainless Steel Alloy for optimal strength and corrosion resistance



**TPSLCOSM**

Multi Tooth
Straight
Stainless Steel
Micro Tissue Plier



**NHDCPVN**

Micro Castroviejo
Curved Diamond
Dusted Needle
Holder



**TP5090**

DeBakey
Perma Sharp™
Tissue Plier



**SHDPV**

Microsurgical
Scalpel Handle

## WHY YOU SHOULD USE AN IMS™ CASSETTE SYSTEM:

By standardizing your instruments with IMS™ System, additional time and revenue is generated while extending the useful life of your HuFriedyGroup Instruments.

Dental professionals state that by implementing IMS™ System, they save 5 to 10 minutes per procedure allowing the staff to spend more time on value-added activities.[3]

In addition, safer transportation of soiled instruments reduces risk of injury and exposure to blood borne pathogens.



**IMN4167**


Large Cassette with Accessory Area, Black

**IM1000**

2-Pack Hinged Instrument Clips

**Learn more about**
**Pink Matters® Surgical Kit**

# ESSENTIAL SURGICAL KIT




**IMN5107**
IMS™ Infinity Series™ **10 Instrument Cassette, Black**
Dimensions: 152 x 184 x 34 mm




**PMNANOSS**
PM Nanos S (Straight),
EverEdge™ 2.0
Perio Knife, Short




**PMNANOAS**
PM Nanos A
(Angulated),
EverEdge™ 2.0
Perio Knife, Short





**SHDPV**
Microsurgical
Scalpel Handle




**TPSLCOSM**
Multi Tooth Straight
Stainless Steel
Micro Tissue Plier




**SPV**
Castro
Microsurgical




**NHDCPVN**
Micro Castroviejo Curved
Diamond Dusted Needle
Holder

1) Pink Matters® is not owned by or affiliated with HuFriedyGroup.
2) When comparing the HuFriedyGroup DLC coated mirror head and handle to the HuFriedyGroup non-coated stainless-steel mirror head and handle.
3) When compared to single instrument decontamination. Based on market survey results. Data on file. Available upon request.

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted.
Pink Matters® is a trademark of The Pink Matters Group.
©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-493/2501

HuFriedyGroup
The Best In Practice